Startup and Shutdown (AeroDR Detector)

Startup

- **1** Start the AeroDR Interface Unit, AeroDR Interface Unit2, AeroDR Generator Interface Unit, AeroDR Generator Interface Unit2, and the access point.
 - · For the startup procedure of the respective devices, refer to "AeroDR SYSTEM/AeroDR SYSTEM 2 Operation Manual".

device.

power.

2 Press the power switch on this

· The startup screen is displayed. · The initialization of the internal system is performed while the startup

· It takes a few minutes to complete the initialization after turning on the

screen is displayed.

· Confirm that each device is turned on.







Patient registration screen





· Upon completion of the initialization, the initial screen (the patient registration screen or the list screen) is displayed, and this device is ready for operation.

3 Confirm the status of the AeroDR Detector.

In case of the wireless connection

If the AeroDR Detector is inserted into the AeroDR Battery Charger/ AeroDR Battery Charger2, remove it.

• In case of the wired connection

Confirm that the AeroDR UF Cable or the AeroDR I/F Cable is connected securely to the wired connection connector of the AeroDR Detector.

AeroDR Detector



4 Start the AeroDR Detector and the AeroDR Battery Charger/ AeroDR Battery Charger2.

· Confirm that each device is turned off.











· Confirm that the AeroDR Detector is ready for use on the screen of this device.





Shutdown

AeroDR Detector



1 Shut down the AeroDR Detector. · Confirm that the AeroDR Detector is turned off.

the access point. Confirm that each device is turned off.

CS-7



DIRECT DIGITIZER



2 Shut down this device.

2-1 Press [System] in the initial screen (the patient registration screen or the list screen).



2-2 Press [Shut down] on the system monitor screen, and press [OK] on the confirmation screen.



• The shutdown screen is displayed, and the shutdown process of this device is carried out. After a while the power is turned off automatically.

3 Turn off the AeroDR Interface Unit, AeroDR Interface Unit2, AeroDR Generator Interface Unit, AeroDR Generator Interface Unit2, AeroDR Battery Charger, AeroDR Battery Charger2, and

* For the operations not described here, please refer to "CS-7 Operation Manual".

Exposure methods (AeroDR Detector)

Exposure

1 Register a patient or select an examination.

- To register a patient, use the patient registration screen.
- To select an examination, use the list screen.

Patient registration screen



- Input patient ID, name, sex, and date of birth.
- · Press [Enter] for each input of items.
- · Press [Perform] after the input is completed.





· Select an examination target to press [Perform].



2-3 Select the Exam. Tag.

- The selected Exam. Tag is added to the Exam. Tags display area. By repeating from step 2-1, Exam. Tags can be added.
- After adding a Exam. Tag, press [Perform].



3 Perform exposure.

· Press the exposure switch of the X-ray device.



set output destination.







- **2** Select the exposure body part.
- 2-1 Select the favorite tab. (Only if there is more than 1 tab)

2-2 Select the exposure body part.



- the image display area.
- The exposed image is displayed in



- When there is more than 1 order, the focus moves for the next order auto-
- matically.



Image adjustment



Using the Image Adjustment tab on the control panel, adjust the image.



- More detailed image adjustment can be performed by displaying the viewer screen.
- The viewer screen can be displayed by press [Viewer] on the exposure screen. Pressing the buttons for [ROI Adjustment], [Parameter], etc. on the control panel image adjustment tab also displays the viewer screen.

Examination completion

When the examination is finished, exposed images are batch output to a pre-

* For the operations not described here, please refer to "CS-7 Operation Manual".

Startup and Shutdown (CR cassette)

Startup



Startup screen



Patient registration screen abcdefg_ pqrstu 456+



• Upon completion of the initialization, the initial screen (the patient registration screen or the list screen) is displayed, and this device is ready for operation.



1 Press the power switch on this device.

- · The startup screen is displayed.
- · The initialization of the internal system is performed while the startup screen is displayed.
- · It takes a few minutes to complete the initialization after turning on the power.

Measures when the REGIUS series do not start automatically?

· Confirm that the power of the REGIUS series is turned on. When it is not turned on, follow the procedures below to turn it on.

1 Turn on the power breaker.



2 Press [operation].



Shutdown



Press [System] in the initial screen 1 (the patient registration screen or the list screen).





[•] The message display area is turned off.

DIRECT DIGITIZER

- **2** Press [Shut down] on the system monitor screen, and press [OK] on the confirmation screen.
- The shutdown screen is displayed, and the shutdown process of this device is carried out. After a while the power is turned off automatically.

When the power of this device is turned off, the connected REGIUS series shifts into sleep.



REGIUS MODEL 110 • The operation lamp blinks slowly.

* Do not turn off the power breaker of the REGIUS series.

* For the operations not described here, please refer to "CS-7 Operation Manual".

Exposure methods (CR cassette/ bar code registration)

Exposure

1 Register a patient or select an examination.

- To register a patient, use the patient registration screen.
- · To select an examination, use the list screen.

Patient registration screen



- Input patient ID, name, sex, and date of birth.
- · Press [Enter] for each input of items.
- · Press [Perform] after the input is completed.







· Select an examination target to press [Perform].

2 Select the exposure body part.

2-1 Select the favorite tab. (Only if there is more than 1 tab)





2-3 Select the Exam. Tag.

- The selected Exam. Tag is added to the Exam. Tags display area. By repeating from step 2-1, Exam. Tags can be added.
- After adding a Exam. Tag, press [Perform].



3 Register a bar code.

· Focusing on the image display area to register CR cassette, read the CR cassette bar code with a bar code reader.







is displayed in the Exam. Tag display area, and the REGIUS becomes ready for reading.

4 Expose and read CR cassette with REGIUS.

- · Press the exposure switch of the X-ray device to expose.
- · After exposing, set the direction of the cassette as shown below to insert it into the slot.

In case of REGIUS MODEL 210





set output destination.

1 Press [**1**] after the exposure completion.





· The exposed image is displayed in the image display area.







DIRECT DIGITIZER

Image adjustment



Using the Image Adjustment tab on the control panel, adjust the image.

- More detailed image adjustment can be performed by displaying the viewer screen.
- The viewer screen can be displayed by press [Viewer] on the exposure screen. Pressing the buttons for [ROI Adjustment], [Parameter], etc. on the control panel image adjustment tab also displays the viewer screen.

Examination completion

When the examination is finished, exposed images are batch output to a pre-

* For the operations not described here, please refer to "CS-7 Operation Manual".

Exposure methods (CR cassette/reading sequence specification)

Exposure

1 Register a patient or select an examination.

- To register a patient, use the patient registration screen.
- · To select an examination, use the list screen.

Patient registration screen



- Input patient ID, name, sex, and date of birth.
- · Press [Enter] for each input of items.
- · Press [Perform] after the input is completed.







(Only	if there	is	more	than	11	tab)

List screen

to press [Perform].

· Select an examination target



2-3 Select the Exam. Tag.

- The selected Exam. Tag is added to the Exam. Tags display area. By repeating from step 2-1, Exam. Tags can be added.
- · After adding a Exam. Tag, press [Perform].



3 Specify the read in sequence.

· Press the Exam. Tag display area in the sequence of exposures to specify the read in sequence.





• The pressed sequence is displayed as [#1], [#2], [#3] and so on in the Exam. Tag display area.

4 Expose and read CR cassette with REGIUS.

-- 🗉 😐

· After exposing, set the direction of the cassette as shown below to insert it into the slot.

In case of REGIUS MODEL 210



· The exposed image is displayed in the image display area.

· CR cassette is read in the sequence of the exposures.

In case of REGIUS MODEL 110









3

set output destination.

1 Press [] after the exposure completion.





















- **2-2** Select the exposure body part.

DIRECT DIGITIZER

Image adjustment



Using the Image Adjustment tab on the control panel, adjust the image.

- More detailed image adjustment can be performed by displaying the viewer screen.
- The viewer screen can be displayed by press [Viewer] on the exposure screen. Pressing the buttons for [ROI Adjustment], [Parameter], etc. on the control panel image adjustment tab also displays the viewer screen.

Examination completion

When the examination is finished, exposed images are batch output to a pre-

* For the operations not described here, please refer to "CS-7 Operation Manual".