


# **Exa® PACS/RIS**

## **Feature Summary**

**Scan and Upload Documents as DICOM Files**

## DICOM Upload

1. **Right-click** a study.
2. Hover over the **DICOM Upload** submenu.



3. Select one of the following options.

<p><b>Upload file</b></p> <ul style="list-style-type: none"> <li>• Upload document (PDF, image) from your computer</li> </ul>	<p><b>Scan to Study</b></p> <ul style="list-style-type: none"> <li>• Scan a document if the study already has images</li> </ul>	<p><b>Scan to Order</b></p> <ul style="list-style-type: none"> <li>• Scan a document if the study doesn't already have images</li> </ul>
---	---	--


## Scan to Study or Order

1. Select the **Doc. Type** (series description).
2. Select **Scan Now**.
3. Select **Accept**.

You can **rotate, flip, and crop** the image. Select **Apply** to save your changes.

## Upload a File

1. Select the **Series Description** (document type).
2. Select **SELECT FILES**.
3. Browse for and select on or more files to upload.
4. Select **UPLOAD**.



5. Confirmation

