

Exa[®] PACS/RIS



New Features and Improvements

34.0.0


©2025 Konica Minolta Healthcare Americas, Inc.

3D Advanced Visualization for MPR


With purchase of a license, 3D Advanced Visualization generates MPR planes and a 3D projection that offer **higher-accuracy rendering and measurements**, and **selectable rendering types** such as Average MIP, minMIP, MIP, parallel, and perspective.



- In any MPR plane:**
- **Rotate and drag** reference lines. All other planes update.
 - Right-click to set the **center point**.
 - Change to Axial, Coronal, or Sagittal plane.
 - Select a **window/level preset**.
 - **Scroll** through images.
 - Use most standard viewer toolbar tools.



In any MPR plane, drag reference lines to adjust the **slab thickness**. Other planes update to match.




In the 3D projection plane:

- Drag the image to **rotate** in 2 dimensions.
- Change to Axial, Coronal, or Sagittal plane.
- Select a **rendering type**, such as Parallel, or Thick and Thin MIP.
- Select an **image type preset**, such as Bond Skin, Lung Film, or CT Fusion MIP.

Viewer > Viewer Settings > Tools > 3D ADVANCED VISUALIZATION

Customize which tools are available when using 3D Advanced Visualization.



View caching status from the worklist


On the worklist, the status of any study caching that is taking place is now indicated with text formatting, and updated in real time.

Dimmed: Uncached *Italics:* Caching in progress **Bold:** Cached


☆ ⓘ	Unread	362	TesA1051	Test, Ashley
☆ ⓘ	Unread	2271840	NE878	Radio, Button
☆ ⓘ	Unread	360	TesA80	Test, Andrew ← Cached
☆ ⓘ	<i>Approved</i>	359	TesB101	Test, Bill ← In progress
☆ ⓘ	Approved	356	TesA1047	Test, Amy
☆ ⓘ	Unread	2239788	RE3434	Re, Mg ← Uncached

Export and import transcription templates

SETUP > OFFICE > Scheduling & Codes > TRANSCRIPTION TEMPLATE




- Export selected or all templates.
- The **Import** dialog now stays open until you close it—import as many templates as you like.




Redesigned transcription template screen

SETUP > OFFICE > Scheduling & Codes > TRANSCRIPTION TEMPLATE




1. **Find** templates quickly with multiple criteria.
2. Identify, search, and sort templates in the list by **Global** availability and **Active** status.
3. **Copy** existing templates to quickly create new ones.



4. **Edit** templates easily.
 - Select **EDIT** to edit the template properties shown.
 - Select [edit icon] to open the main template editing screen to modify contents and layout.

Dynamic selection of markets, facilities, modalities, and modality rooms

In various parts of the program, new controls are now available to search for and multi-select markets, active facilities, modalities, and modality rooms.



- Selecting **markets** auto-selects corresponding facilities in the **Facilities** list, and vice versa.
- Clear search/filtering criteria
- Search/filter within the list
- Select all
- Items appear in alphabetical order
- Selecting **modality rooms** (not shown) auto-selects corresponding items in the **Modalities** and **Facilities** lists, and vice versa.


Add Patient Guarantors from the Edit Study screen

Edit Study > PATIENT GUARANTOR > ADD



Remove OpalRAD Annotations from teaching studies


Worklist > right-click a DICOM study > Images > Create Teaching Study



Send scanned documents as images or PDFs


SETUP > OFFICE > DICOM > AE TITLE > ADD

You can configure AE's to send scanned documents as DICOM images or encapsulated PDFs.




NPI added to the provider pop-up window on the worklist

Now you can view the NPI of the provider in the worklist pop-up window.



More control over passwords

SETUP > OFFICE > COMPANY > PASSWORD MGT.



Require (or not require) passwords to contain specified numbers of the following:


- **Min./max. characters**
- **Uppercase & lowercase letters**
- **Numbers**
- **Symbols & special characters**

Users Cannot Reuse the Last [] Passwords

Specify how many previous passwords users may not re-use.


Guidance when changing your password

Burger  > Globe  > CHANGE PASSWORD



Reorder charges on a claim

On the claim editing screen, you can reorder charges to control how they appear on the finished claim. Use the  and  buttons to move charges up and down in the list.



Edit : Autotest_Lastaiytczqwo, Autotest_Firsti62gw4opg

CHARGES CLAIMS INSURANCE ADDITIONAL INFO BILLING

PREV

				NDC	Date	Accession No.	CPT
	+	×			06/29/2024	1605	112
		+	×		06/29/2024	1605	768
	+	×			06/29/2024	1605	768

NEW CHARGE

Sort patient claims by policy and group number

Billing > Claims > Claim Inquiry > Patient Claims

When you want to group personal injury cases or episodes together and then view or print the selection on the Patient Activity Statement, you can sort patient claims by policy number and group number.



Patient Claims: Test, Joe

PREVIOUS NEXT

Billing Provider: All

Unapplied Payments: \$0.00

All Activity Activity by Date


Include Billing Provider Info
 Include Third Party Payers
 Include All Comments

PRINT FAX

CLAIM NO.	STUDY DATE	BILLING FEE	TOTAL ADJUS	TOTAL INSURANCE PA	TOTAL PATIENT PAYM	BALANCE	CLAIM STATUS	CURRENT RESPONSIBILITY ASSIGNED	POLICY NUMB	GROUP NUM
<input type="checkbox"/> 1028	10/02/2024	\$0.00	\$0.00	\$70.00	\$0.00	(\$70.00)	Pending Payment	AARP MEDICARE COMPLETE	18383	96666
<input type="checkbox"/> 1027	10/02/2024	\$50.00	\$0.00	\$50.00	\$0.00	\$0.00	Paid in Full	NATIONWIDE HEALTH PLANS	4477474	1549

New merge fields for study forms

You can now add merge fields to study forms that populate with information from all studies in an order. For multi-study orders, patients using Patient Kiosk no longer need to fill out a separate form for each exam.



New Merge Fields:

- Ref_PhyFullName_order
- TechnologistFullName_order
- AttorneyFullName_order
- Ord_PhyFullName_order
- AccessionNo_order
- Modality_order
- ModalityRoom_order
- StudyDate_order
- StudyDescription_order
- StudyFlag_order
- StudyReceivedDate_order
- ExamPrepInstructions_order
- AppointmentTypeCode_order
- AppointmentTypeDescription_order
- StudyICDCode_order
- StudyICDDescription_order

Add modifiers to charges in the Technologist screen

Technologists screen > Add Charges


You can now add modifiers directly in the Technologist screen. This is useful when claiming reimbursement for unused drugs or contrast as part of NDC claim submissions.



Configure document types for AEs

SETUP > OFFICE > DICOM > AE TITLE > ADD

Select which documents to DICOM-send on each AE.




1. Clear the **Ignore Scan Docs (Send)** checkbox.
2. In the **Document Types** dropdown list, select which documents to include in any DICOM-send operation.

Alberta claims: Select from all active referring provider contacts

Billing > Claims > [Edit claim] > Referring Provider

When adding a referring provider to a claim, you can now select from all active contacts (excluding system providers) rather than just the primary contact.



Alberta: Enter the WCB billing number for radiology providers

SETUP > OFFICE > Providers & Resources > RESOURCE > ADD

To meet new WCB requirements, you can enter a WCB billing number for your radiology provider resources. The number is automatically included when submitting claims.

The screenshot shows a form for adding a new resource. The 'Type' dropdown is set to 'Provider-Radiology'. The 'WCB Billing Number' field at the bottom is highlighted with an orange box. Other fields include Code, Name (First Name, MI, Last Name, Suffix), Title, Taxonomy Code, Dr. Office Name, Federal Tax ID, EIN No., ETIN No., Medicare UPIN, and an 'Inactive' checkbox.

Dispatching Dashboard: Reset all studies in order


DISPATCHING DASHBOARD > ASSIGNED

On the ASSIGNED tab, you can select **Reset All Studies in Order** from the shortcut menu to reset studies to the Scheduled status (requires the Reset All Studies in Order user right). Not available for orders in a transcription-related status.

The screenshot shows the Dispatching Dashboard with the 'ASSIGNED' tab selected. A table lists orders with columns for 'CREATED DATE', 'SCHEDULE DATE', and 'REQUESTING DATE'. A context menu is open over a row, with the 'Reset all studies in order' option highlighted by an orange box. Other menu options include 'Edit Order', 'Patient Information', 'Change Status', 'Stat', and 'Priority'.

Dispatching Dashboard: Sort by dispatch number


In Dispatching Dashboard you can now sort orders by the **first three digits** of the dispatch number.



Dispatching Dashboard: Add to Open Dispatch

DISPATCHING DASHBOARD > TO BE DISPATCHED > [Right-Click] > Dispatch.


Adding orders to existing dispatches is now easier in the Dispatch dialog. Also, now you can add orders to open dispatches whether or not they are associated with ordering facilities.



1. In the upper-left dropdown list, select **Technologist** or **Vehicle**.
2. In the **Add to Open Dispatch** dropdown list, select the existing dispatch to add to.
3. Modify other settings if needed, and select **DISPATCH**.

Improvements to the Schedule Rules screen


SETUP > OFFICE > Scheduling & Codes > SCHEDULE TEMPLATES



1. Search for schedule rules by **repeating rule date** or **date/time rule date**.
2. Sort rules by whether they appear as blocks on the schedule book.
3. Open a summary screen to view rule details.

More information in Mobile Rad

Select and view more information in Mobile Rad, including **patient alerts**, **ICDs**, **study notes**, and **reasons for study**.



Manually set eligibility

Worklist > [Shortcut menu] > Actions > Manually Set Eligibility

From the worklist, you can now manually set eligibility of a study to **Active** or **Inactive**, or select to **Clear** the manually-set status.

When **Active**, a check mark appears in the **ELIGIBILITY** column.

View and edit NDCs in claims

Billing > Claims > [Double-click a claim]

When editing charges on a claim, you can now view and edit NDCs, if available.

			NDC	Date	Accession No.	CPT Code	CPT D
	⊕	⊗	0270-5164	11/15/2024	19118	77076	RADI
	⊕	⊗	+	11/15/2024	19118	Select	Sele

NEW CHARGE

View approved reports in Addendum-related statuses

You can now view approved reports of studies in **Draft – Addendum** and **Addendum – Transcribed** statuses.

+ [Icons]	Unread	11/20/2024 1:30 AM EST	Lcspa
+ [Icons]	Draft - Addendum	11/20/2024 1:30 AM EST	Ulspat
+ [Icons]	Unread	11/20/2024 1:00 AM EST	Lcspa

Simplified navigation in Patient Portal

Go straight to where you need with the redesigned navigation pane.

The screenshot shows the 'PATIENT INSURANCE' section of a patient portal. On the left, a navigation pane is highlighted with an orange border and contains the following items: Appointments, Request Appointment, Calendar, Insurance, Guarantors, Authorized Representatives, Account, and Activity Log. An orange arrow points from the 'Request Appointment' item to the main content area. The main content area displays insurance details for 'AARP MEDICARE COMPLETE' with a policy number '1234'. The text 'Simplified, centralized navigation pane' is written in orange next to the arrow.

New user right for assigning study filters

If granted the new **Assign Study Filters** user right, non-administrators can assign worklist study filters to other users.

Study Filter

Filter Name * Show Encounters Only

Joined Filters

Filter Order * Is Private

Display as a Tab Display in Dropdown Show Only Exceeding Max Time

Show DICOM Studies Only Show RIS Orders Only Show Assigned Studies Only

Deleted Show Pre-Orders Only

Assign

Configure Philips Speechmikes

In **Exa Trans**, on the **Options** tab, you can map functions to **Philips Speechmike** buttons.

ExaTrans

Transcription Documents Prior Reports Monitor **Options**

Application
Theme

Show Hypothesis
 Pin Hypothesis
 Show "Enable Approvals" Checkbox

Editor
 Keep Delimiters
 Replace Entire Word

Microphone
 Hold Button To Record

Button Mappings

Button	Action	Parameter
FastForward	Previous Field	
Rewind	Next Field	
Command	Approve and ...	
Eol	Approve	
Insert	Pre-Approve	
Instruction	None	
Function1	None	
Function2	None	
Function3	None	
Function4	None	

Add modality room notes

SETUP > OFFICE > Office > MODALITY ROOM > MODALITY ROOM NOTES

You can now add notes for each of your modality rooms that are available to view on the schedule book.

MODALITY ROOM MODALITY ROOM NOTES

SAVE SAVE & CLOSE

Modality Rooms * ANGIE CT XA 🔍

Start Date * 01/01/2024 📅

End Date 01/31/2025 📅

Notes * Preferred modality room for CT XA.

1. Select a **modality room**.
2. Set **Start** and **End** dates.
3. Type the **notes**.

☰ Schedule Book

RADIOLOGY FACILITY : ANGIE FACILITY MYFACILITY FILTER

Angie CT XA (00ang) Angie MR (00ang)

07:00 am

07:15 am

07:30 am

07:45 am

Notes


Preferred modality room for CT XA. ✎

On the **schedule book**, if the blue information (i) button appears, you can select it to view the notes.

Clear critical findings

Worklist > [shortcut menu] > Actions > Critical Findings



You can now clear any previously added critical findings from a study.



Manage third-party apps

SETUP > OFFICE > User Management > USERS > [open a user record]


Now administrators can easily manage a user's third-party apps directly in their user record.



Assign dates to report templates


SETUP > OFFICE > Scheduling & Codes > REPORT TEMPLATE > ADD

You can now assign a start and end date to your report templates. This makes the template available only to studies whose study dates fall within that date range.



Display of number of reasons per category


SETUP > OFFICE > Office > REASON CODES



Inactivate patient guarantors, and manage from Edit Study

PATIENT > PATIENT GUARANTOR

Now you can manage patient guarantors in both the patient chart and the Edit Study screen. Only one guarantor can be active at a time.



Add primary and additional technologists to individual studies

You can assign the primary technologist to individual studies (rather than orders) directly in the **Edit Study** and **Technologist** screens. Plus, you can now assign up to 10 **additional technologists** to be included in the Technologist Productivity operational report.

