



KONICA MINOLTA

ImagePilot Sigma

with TCO

A Total Cost of Ownership CR Solution

- Simplify the transition to digital
- Manage long term cost

For practitioners converting from film to digital or upgrading their CR systems, Konica Minolta's ImagePilot Sigma with Total Cost of Ownership (TCO) coverage is an innovative solution combining state-of-the-art CR technology and long term support services into one simple and cost effective purchase.

- Outstanding Image Quality
- Versatile Workflow
- Total Cost of Ownership Management
- Help Secure Patient Data



Outstanding Image Quality

ImagePilot Sigma achieves outstanding image quality consistently and easily with AutoPilot automatic image processing software. This simple, single click technology means there is no need to specify body part, exam, or view. AutoPilot displays the optimized diagnostic view by compensating for technique variability and by learning user preferences. The result: incredibly consistent image quality for confident diagnosis and satisfied users.

Versatile Workflow

The ImagePilot Sigma is an All-in-One solution that provides simple image acquisition with PACS expandability. With the same intuitive user interface at each workstation, ImagePilot Sigma enables image review or CD output anywhere in the office. Easily access advanced measurement and annotation tools from the X-ray, examination, or consultation room. And with extensive connectivity, ImagePilot will integrate with most EMR software and PACS systems to help simplify office workflow and improve staff efficiency.

Giving Shape to Ideas

Total Cost of Ownership Management

ImagePilot Sigma TCO extends customer satisfaction and enhances economic value with long term coverage that protects your investment and minimizes ongoing costs. Our all-inclusive TCO plans cover remote technical support, replacement parts, software upgrades, cassette repairs, and Informity, so you can stay focused on your patients and not worry about technology.

Help Secure Patient Data with Informity

Help satisfy HIPAA requirements by securing your protected health information with Informity, Konica Minolta's cloud-based subscription service. Informity automatically backs up patient and system data to a secure datacenter, mitigating the risk of lost patient data and providing seamless recovery in the event of a disaster.

Hardware Components

Imaging Station	ImagePilot Sigma ACQ	ImagePilot Sigma 25, 45, 60	ImagePilot Sigma Mobile
Server Platform	Desktop	Desktop	Laptop
Hard Drive Capacity	500GB (1x500GB)	500GB (1x500GB)	320GB (1x320GB)
RAM	4GB	4GB	4GB
Display Options	19" LCD	19" LCD or 23" Eizo Multi-Touch	15.6" Built-in
UPS for Imaging Station	Optional	Optional	Optional
REGIUS CR Reader	REGIUS Sigma II	REGIUS Sigma II	REGIUS Sigma II
Throughput (14x17 cassettes per hour)	45 or 60	25, 45, or 60	45 or 60
Backup/Extended Storage			
DVD Backup Tool-For Disaster Recovery	Optional	Yes	Yes
Backup NAS, 4TB or 8TB RAID 5	Optional	Optional	Optional
Extended Storage NAS, 4TB or 8TB RAID 5	Optional	Optional	Optional

Total Cost of Ownership Coverage

TCO Benefit	ImagePilot Sigma ACQ	ImagePilot Sigma 25, 45, 60	ImagePilot Sigma MOB
Three years of Parts and Technical Support	Yes	Yes	Yes
Three years of Cassette Repair one per year, \$150 copay	Yes	Yes	Yes
Three years of Software Upgrades	Yes	Yes	Yes
One Year of Informity Data Bank (annual allotment)	No	Informity 1,500 65 GB ~1,500 exams	Informity 1,500 65 GB ~1,500 exams
One Year of Informity Collaboration Space	No	100 MB About 5-10 exams at any time	100 MB About 5-10 exams at any time

As the MD Buyline #1 rated CR manufacturer since 2005, Konica Minolta delivers performance, reliability, and durability with its CR solutions. Benefit from quality solutions and service from a dependable company and a valuable partner.

Konica Minolta's Commitment to Customer Satisfaction

At Konica Minolta, your complete satisfaction is our number one goal. From pre-installation through the lifetime of your investment, our professional Reseller Network provides the technical expertise you need. Project management, Applications Training, and Field Services are different teams of Konica Minolta Trained Reseller professionals that work together to maintain your continued satisfaction over the lifetime of your investment.

Maximize your investment in CR equipment by selecting the right Customer Satisfaction Agreement to ensure the ongoing operational excellence and productivity that your department needs.



KONICA MINOLTA

Konica Minolta Medical Imaging

411 Newark Pompton Turnpike
Wayne, NJ 07470

Tel: (973) 633-1500 Fax: (973) 523-7408

Website: konicaminolta.com/medicalusa