
## Imaging analytics for better decisions, sooner.






# Maximize system utilization. Minimize service interruptions. Optimize personnel performance.

AeroRemote<sup>®</sup> Insights is a unique monitoring and analytics tool that makes it far simpler and more efficient for you to manage digital radiography assets, run a more productive imaging department and deliver a better experience for patients.

A cloud-based subscription service available with Konica Minolta Digital Radiography Systems, AeroRemote Insights gives you up-to-date metrics on procedure volumes, staff performance, and system health, facilitating quick, confident decision making on critical aspects of your department's performance.





#### **Track Departmental Productivity**

Gain deeper understanding of how, when and by whom systems are being used—at a glance.



#### **Evaluate User Performance**

Assess user accuracy, productivity, and exam dose patterns to help evaluate performance and identify training opportunities.



#### **Maintain Optimal Performance**

Stay on top of system performance, maintenance requirements and service contracts, with a click.




#### **Respond to Critical Events**

Get immediate notification about critical events, such as panel drops, so you can respond quickly and minimize service interruptions.

### Interactive Dashboard Brings Analytics to Life

AeroRemote<sup>®</sup> Insights collects valuable data daily from your DR systems<sup>1</sup> — from exam volumes to rejected image reasons and exam exposure levels. This data is automatically integrated into intuitive presentations that help you streamline department management and decision-making.




#### **Productivity**

Learn at a glance how each hospital, department, system or user is performing—by time or exam type. Gain insight into asset use and staff training opportunities.

Operator Username	Exposure	Average S Value	Average El Value	Average Di Val
CRR	7,177	473	298.33	-1.
DC	678	496	287.61	-1.
DCJ	6,489	485	253.16	-1.
DMM	1,176	485	269.42	-1.
KAC	618	514	328.31	-12
KLC	5,615	502	324.48	-3.
SMT	2,069	237	657.57	0.
TJW	10,951	489	232.03	-1.7
VLM	9,150	456	296.50	-3.
VRW	3,418	413	357.73	-0.


#### **Dose Metrics**

Monitor exposure indicators by user, system or exam view. Identify where techniques can be adjusted to optimize dose.




#### **Aggregate Comparisons**

Compare your facility's performance to the AeroRemote Insights community average, including reject rates.




#### **Staff Performance**

Track image rejects to identify causes of repeated exposures. Identify need for training or system maintenance.



#### **System Health**

Gain instant access to service contract status, system events and AeroDR panel calibration schedules. Minimize downtime by addressing issues before they become problems.



#### **Panel Handling Tracking**

Assess AeroDR panel handling by time of day, location and user. Identify where to take action to reduce the financial risk of accidental damage and downtime.

## Unlock value department-wide.









Manage your DR business from anywhere - AeroRemote Insights analytics dashboard and monitoring service is available for most AeroDR models and is compatible with many iOS and Android mobile devices.

#### Easy to use and feature-rich:

- · Easily export charts and data to spreadsheets or PowerPoint
- Email presentations to external colleagues
- Drill down to additional data or views with a simple click
- Create your own filters for viewing data by a specific criteria
- · No patient specific data is ever transmitted
- Data updates daily

#### View data by:\*

- Facility name
- · Room/location
- Time
- Username
- Anatomical view
- \* filtering may vary depending on the metric or chart being viewed

#### **Available metrics:**

- Exam and exposure counts
- Rejected image statistics (counts, rates, reasons)
- Exposure indices
- System service contract status
- Operating system warnings
- Impact (drop) counts

#### **System requirements:**

- CS-7 Version 1.35 or later; Ultra software version 4.32 or later
- Konica Minolta R3 connection agent for each CS-7 or Ultra station
- Internet access for each CS-7 or Ultra station
- Use unique logins for maximum utility of the usage screens



AeroRemote Insights is included in Konica Minolta's comprehensive Blue Moon Select service agreements or is available as a stand-alone annual subscription.



Watch the AeroRemote Insights animated video to learn more. konicaminolta.com/medicalusa/insights

1. Compatible with CS-7 or Ultra acquisition software (may require software upgrade). Not all panel models detect drops.

© 2020 Konica Minolta Healthcare Americas, Inc.





konicaminolta.com/medicalusa