**AeroDR<sup>®</sup> GL-P** is a high dynamic range glassless digital detector solution for general radiography.



BEND RESISTANCE Up to 286 lbs.



WATERPROOF IP56 Compliant

**STRONGER** Holds up to 881 lbs.



60% Less Grip Strength Needed



## The AeroDR<sup>®</sup> GL-P advantage.

The simple way to upgrade to DR. Our affordable and simple-toinstall retrofit solutions that add the benefits of Digital Radiography capabilities to your existing x-ray.

Konica Minolta delivers solutions that provide clinical and economic value to meet your needs. The AeroDR Glassless Flat Panel Detectors are our lightest, strongest and most sensible solution to help you improve efficiency and reduce user fatigue without compromise.

## Increased sensitivity by replacing the glass substrate with a thin film transistor (TFT)

AeroDR Glassless Flat Panel Detectors use a thin film instead of a glass substrate to achieve lighter weight. Removing the glass reduces the distance between the scintillator and the TFT to increase sensitivity. The electronics are optimized to reduce electrical noise. Additionally, the AeroDR GL-P provides a High Dynamic Range option that provides a 200µ resolution for a wider range of tone variations to visualize tissue.

## Superb image quality - Glassless standard resolution imaging

- High-dynamic range imaging at 200µm resolution
- Power cell life up to 4.3 hours on full charge (1 LiC)
- 14" x 17" panel weighs 3.9 lbs/1.8 kg super lightweight, including the capacitor
- IPX56 dust and liquid resistance (with integrated power supply)
- 4.0 sec AED with AeroSync®
- Image area 348.8 x 425.6mm
- External dimensions 15.1"x18.11"x 0.6"
- Matrix size 1,744 x 2,128 (200µ)
- Image preview in ~2 seconds
- Diagnostic view in 4 -7 seconds

High dynamic range takes advantage of REALISM image processing to enhance tone variation to visualize soft tissue with minimal impact to resolution at normal viewing ratios.







For more information about these products, please contact your Konica Minolta Sales Representative.

© 2023 Konica Minolta Healthcare Americas, Inc.



Konica Minolta Healthcare Americas, Inc. 411 Newark Pompton Turnpike Wayne, New Jersey 07470 Tel: +1 (973) 633-1500 km.marketing@konicaminolta.com healthcare.konicaminolta.us

M2090 1222 RevB