



Exa® PACS/RIS

1.4.31_P2

Customer Release Notes

© 2022 Konica Minolta Healthcare Americas, Inc.

Table of Contents

New Features	3
Customer Requested Improvements	3
Other Improvements	3

New Features

Summary

- No new features in this release.

Customer Requested Improvements

Summary

- In the Alberta configuration, when editing claims, now only the primary contact of a referring provider can be selected.
- Users can now deactivate insurance policies, and easily differentiate between active and inactive policies in various parts of the program.
- Further improved the way insurance policies are made active or inactive, how the Valid To date is determined, and how that information is communicated to the user.
- Improved the way Exa PACS/RIS handles study statuses of studies being editing by multiple simultaneous users.
- Diagnosis codes are now sent to third party billing systems via HL7 in the correct order.

Other Improvements

Summary

- When constructing URLs to launch third-party applications, additional variables are now available for patient name, DOB, sex, and study description.
- When using MPR, the program now notifies the user if a high gantry tilt causes measurements to become approximate.
- Various improvements to performance and the UI.