

exa-PACS • exa-RIS

Exa® PACS/RIS

1.4.31_P3

Customer Release Notes

© 2022 Konica Minolta Healthcare Americas, Inc.

Table of Contents

New Features	3
Customer Requested Improvements	3
Other Improvements	3

New Features

Summary

- No new features in this release.

Customer Requested Improvements

Summary

- The Exa platform now generates Opal Forwarder licenses for active physical adapters, and only requires such licenses to run.
- PR (GSPS) DICOM objects sent from Intelrad are now recognized by the Exa PACS/RIS viewer as DICOM overlays that can be turned ON and OFF.
- Improved the way the Exa PACS/RIS viewer parses DICOM tags to determine which annotations to display.
- When changing a study's appointment types, users can now copy selected authorization information from old to new CPTs.

Other Improvements

Summary

- When a user selects a 2D image annotation to jump to the corresponding 3D instance, Exa PACS/RIS now displays the corresponding annotations on the 3D instance.
- Users can now configure a color in which to display 3D CAD SR markings when hovering over them.
- Various improvements to performance and the UI.