

exa-PACS • exa-RIS

Exa® PACS/RIS

1.4.31_P4

Customer Release Notes

© 2022 Konica Minolta Healthcare Americas, Inc.

Table of Contents

New Features	3
Customer Requested Improvements	3
Other Improvements	3

New Features

Summary

- Users can now configure whether to send faxes upon receipt of inbound HL7 ORU messages.
- Users can now configure whether to output HL7 messages, faxes, emails, and printouts based on HL7 global, study status, and routing rule triggers.

Customer Requested Improvements

Summary

- Exa Dictation can now be used with any version of Exa PACS/RIS.
- When scheduling appointments back to back, users can now specify a maximum delay between appointments of Same Day, which returns all available timeslots on that day with 0 or more hours of delay.
- Schedule rules names can now include up to 511 characters.
- Added a workflow for creating quick blocks from the schedule book for users with both quick block and schedule block rights.

Other Improvements

Summary

- The DICOM and Opal Listener services now provide compression and transcoding functionality for received images.
- Opal Import now supports importing of series in which DICOM tag 0028;0004 is set to MONOCHROME2.
- Various improvements to performance and the UI.