

exa-PACS • exa-RIS

# **Exa® PACS/RIS**

## **1.4.32\_P1**

### **Customer Release Notes**

© 2022 Konica Minolta Healthcare Americas, Inc.

## *Table of Contents*

New Features .....	3
Customer Requested Improvements .....	3
Other Improvements .....	3

## New Features

### Summary

- Users can now add and delete settings for encrypting study URLs, and can generate corresponding encryption keys.
- The Exa PACS/RIS viewer can now be opened by using an encrypted URL from an EMR or other external source.

## Customer Requested Improvements

### Summary

- When entering an adjustment payment and changing the payer, Exa PACS/RIS no longer erroneously changes the payer to the patient.
- Users can now set a span sensitivity of Default, which was the 0/Auto sensitivity used in version 1.4.29.
- Users can now search for patients by middle initial.
- In the HL7 Sender Log screen, the Retry All button was removed.

## Other Improvements

### Summary

- Various improvements to performance and the UI.