

## Quality Tools Integrated with Radiologist Workflow.

## Exa® Peer Review

Exa Peer Review facilitates the feedback process from assignment through assessment and reporting. Peer Review is an add-on module integrated into the radiologist's workflow.\* Tools are provided to effectively review, analyze and take action to improve performance. Patient data is anonymized to eliminate biases and provide objective feedback.


Peer Review is designed for a personalized experience for radiologists. Robust flexibility is provided in the assignment of studies that are presented as a tab in their worklist.


Assignments can be based on a combination of:

- Percentage of yearly approvals
- Facilities
- Modalities
- CPT codes
- · Appointment types

Studies can also be assigned manually.

Double-clicking on an assigned study launches the study and a multi-paneled review screen in which to assign a level of agreement, provide comments and review prior reports.





\*Exa Peer Review requires an additional fee.



Konica Minolta Healthcare Americas, Inc.

2217 U.S. Highway 70 East Garner, NC 27529 Phone: 800.366.5343

konicaminolta.com/medicalusa

M1599 0220 RevA