

exa-PACS • exa-RIS

Exa® PACS/RIS

1.4.32_P4

Customer Release Notes

© 2023 Konica Minolta Healthcare Americas, Inc.

Table of Contents

New Features	3
Customer Requested Improvements	3
Other Improvements	4

New Features

Summary

- As a general rule, when reports are linked, activities related to retaining the STAT level or recording the TAT of a main study now also apply to any of its linked studies.
- Exa PACS/RIS now supports integration with Visage RIS.

Customer Requested Improvements

Summary

- For billing rules of Exam Type Override and Exam Type, percentage-based fee overrides in billing rules are now based on the billable fee schedule assigned to the ordering facility if one exists, and the default fee schedule otherwise.
- In the Alberta and Manitoba configurations, users can now manage fee schedule groups.
- In the Alberta and Manitoba configurations, users can now assign fee schedules to fee schedule groups.
- In the Alberta and Manitoba configurations, users can now create multiple Default and Self-Pay fee schedules.
- In the Alberta and Manitoba configurations, when assigning fee schedules to insurance providers, users can now filter available fee schedules by fee schedule group.
- In the Alberta and Manitoba configurations, when assigning fee schedules to an ordering facility, users can now filter available fee schedules by fee schedule group.
- In the Alberta and Manitoba configurations, when assigning fee schedules to a referring provider, users can now filter available fee schedules by fee schedule group.
- In the Alberta and Manitoba configurations, when assigning fee schedules to a facility, users can now filter available fee schedules by fee schedule group.
- In the Find Slots screen, while in the process of scheduling appointments back to back, users can now restart the process to edit the appointments in selected timeslots without having to start over.
- KM staff can enable the Fee Schedules by Group and Dates feature, which allows users to choose multiple fee schedules from a group.
- Modifiers are no longer automatically applied to setup charges for multiple patient seen, and setup fees are no longer split between multiple patients.
- The configuration option to send approved reports by Fax to an ordering facility is now available only when editing the ordering facility's primary contact/location.
- When configuring ordering facilities, the following information is now managed at the contact/location rather than the main OF: phone number, fax number, email address, and report password.
- When rescheduling an exam and copying its authorization, the expiration date now remains unchanged.

Other Improvements

Summary

- For CD burning, CD/DVD-R and CD/DVD+R type discs are now supported.
- In the Alberta configuration, Exa PACS/RIS now considers the From and To dates of the fee schedule to determine billing fees.
- Logs now record what type of compression is used when burning a CD.
- Made improvements to better handle large volumes of continuously inflowing data, especially to avoid duplication of patient charts and studies.
- OpalPrefetchListener is now called Exa Burner.
- When burning a disc, Exa Burner shows the available space, percentage of space to be taken by the burned studies, and the amount of buffer space used during burning.
- When linking studies with different transcription templates and approving the main study, both main and linked studies now use the template of the main study.
- Various improvements to performance and the UI.