

exa-PACS • exa-RIS

# Exa® PACS/RIS

## 1.4.32\_P5

### Customer Release Notes

© 2023 Konica Minolta Healthcare Americas, Inc.

## *Table of Contents*

New Features .....	3
Customer Requested Improvements .....	3
Other Improvements .....	3

## New Features

### Summary

- Directly from the worklist, users can now create new orders and studies while opening Ultra Acquire to add images.
- In Symmetry PACS, KM staff can now show or hide a button on the worklist that opens Acquire, and buttons for new studies and adding images.
- In Symmetry PACS, users can now open Ultra Acquire from the worklist, or open it specifically in Acquire mode to add images to the current study.

## Customer Requested Improvements

### Summary

- For verifying eligibility from Imagine, Exa PACS/RIS now supports duplicate Imagine provider IDs.
- In OpallImport, AE Titles are now only available to select if the user has access to them and their associated facilities.
- In the Alberta configuration, users can now configure a Query/Retrieve AE to search for patients by ULI or Exa PACS/RIS account number, and to search for all patients or for Exa PACS/RIS patients only.
- In the Alberta configuration, users can now perform query/retrieve of the AHS database by ULI or Exa PACS/RIS account number, and can select to return all patients or only patients in Exa PACS/RIS.
- In the Alberta configuration, when claim charge dates are modified, the fee schedule that matches the updated date is now used to determine fees.

## Other Improvements

### Summary

- Various improvements to performance and the UI.