Universal DR HEALTHCARE

Digital Radiology, Simplified.





ULTRA Image Acquisition Software

The Universal DR system is an economical solution that simplifies the transition from analog to digital

- The combination of advanced analog-to-digital imaging and handling features simplifies the operation and use for radiography in general
- The wireless CsI detector, in conjunction with the imaging processing software ULTRA, delivers exceptional image quality and the right combination of advanced features in simplified software for all general radiography applications
- It features Automatic Exposure Detection (AED), enabling quick installation and operational readiness process, without requiring adjustments or connections to the X-ray equipment
- The Universal DR wireless detector is available in 14" x 17" and 17" x 17" cassettes and is compatible with any standard X-ray system



ULTRA Image Acquisition Software

The combination of the Universal DR wireless detector with customizable ULTRA image processing software makes for an ideal and straightforward DR upgrade or retrofit DR solution. The intuitive user interface, which includes an intelligent workflow design for managing exams on a single screen, enhances the overall image workflow and saves time for both the user and the patient.

Universal DR simplifies the transition from analog to digital radiography and is suitable for both new X-ray systems and upgrading existing X-ray rooms.

Main features of the ULTRA software

- Connectivity compatible with the DICOM standard and RIS/PACS
- Procedure code mapping tool
- Study annexes
- · Free text annotation
- · Automatic masking
- · Study list filter
- History
- · Image zoom
- Grid suppression
- Enabling HIPAA compliance (password protection, screen lock)
- Foldering (to merge and move studies)
- DICOM Store
- DICOM modality work list (MWL)

AeroRemote® Insights


The cloud-based AeroRemote Insights monitoring service is an easy-to-use analytical tool that assists institutions in effectively managing their digital radiography assets, and is available with Universal DR systems.

Access real-time metrics on procedure volumes, staff performance and system health to make quick and informed decisions. Visualize daily data with user-friendly graphs for assessing productivity, identifying training opportunities and optimizing asset management.

You can select AeroRemote Insights with Blue Moon Select service contracts or purchase it separately as a subscription.







BLUE MOON LIFECYLE PLANS

Blue Moon Lifecycle Plans always provide you with the necessary technical support, helping you avoid unexpected expenses. We consistently introduce new features to our plans to offer even more options. Minimize downtime, maximize productivity and eliminate risk with service solutions that meet your needs.



© 2023 Konica Minolta Healthcare Americas, Inc



Konica Minolta Healthcare Americas, Inc.

411 Newark Pompton Turnpike Wayne, New Jersey 07470 Tel: +1 (973) 633-1500 km.marketing@konicaminolta.com

healthcare.konicaminolta.us

M1528 1023 RevC