

AeroDR® PDR Digital Detector Solution for General Radiography.
Value You Can Trust. Performance That Lasts.



BUILT-IN POWER

4.1 HOURS, NO REPLACEMENT



SHARP RESOLUTION

REVEAL DETAILED IMAGES



QUICK CHARGE

15 MINUTE FAST BOOST



SHARPER IMAGES

48% DQE FOR CLEAR,
HIGH-QUALITY IMAGES



LIGHTWEIGHT

LESS THAN 5.5 LBS



CLEANER, SAFER USE

ANTIBACTERIAL SURFACE



FAST WORKFLOW

1S PREVIEW, 4S CYCLE



Built Tough



BEND RESISTANCE

Handles 287 lbs bending force



DROP RESISTANT

Tough at 1.2 m drops



DUST & WATER RESISTANCE

IP56 for rough use



HIGH LOAD

Supports 882 lbs

The AeroDR® PDR. Reliable. Affordable.

AeroDR PDR Wireless Detector delivers high-performance imaging and consistent image quality through advanced processing. As a cost-effective member of the AeroDR family, it balances affordability and performance to meet both clinical and operational needs.

Upgrade Made Easy. A smart path to digital radiography with a versatile, value-focused solution that integrates easily with existing X-ray systems. It helps support clinical confidence while addressing workflow needs and budget considerations.

Next-Generation Features Delivering Value and Quality

Designed for today's clinical demands, this next-generation detector combines smarter electronics, advanced materials, and a refined design for stronger overall performance. Upgraded circuitry reduces image noise and improves consistency, while a more efficient scintillator layer enhances light conversion, producing clearer images at a lower dose. A reinforced aluminum seal increases durability for everyday use. Together, these improvements deliver sharp, reliable images with less downtime, making it a dependable, value-driven solution for facilities that demand performance without compromise.



AeroDR PDR 14"x17" - Key Benefits at a Glance

Workflow Efficiency

- Fast image delivery supports smoother exams and less patient waiting
- Onboard power keeps downtime to a minimum — no battery replacements needed

Design & Durability

- Lightweight and rugged — easy to use, tough enough for busy clinical environments
- Sealed against dust and water, with an antibacterial surface for added hygiene

Imaging Performance

- High-resolution imaging with advanced processing — see fine details clearly, even at lower dose
- Delivers consistent, reliable results that support confident clinical decisions

Trusted Support, Reliable Partnership

Backed by extensive experience and dedicated service, our team provides responsive support and long-term reliability—helping healthcare facilities focus on delivering quality patient care with confidence.



For more information about these products, please contact your Konica Minolta Sales Representative.

© 2025 Konica Minolta Healthcare Americas, Inc.



KONICA MINOLTA

Konica Minolta Healthcare Americas, Inc.
411 Newark Pompton Turnpike
Wayne, New Jersey 07470
Tel: +1 (973) 633-1500
km.marketing@konicaminolta.com
healthcare.konicaminolta.us

M2373 0725 RevA