

DIGITAL RADIOLOGY. YOUR WAY.



Imaging software supports detailed bone and soft tissue visualization.

THE NEXT EVOLUTION IN DIGITAL RADIOGRAPHY

The KDR® FMT and KDR FMT+ represent the next evolution in floor-mounted digital radiography from Konica Minolta. Designed to deliver exceptional image quality, reliability, and workflow efficiency, these systems build on the trusted performance of the KDR product line—making them the preferred choice for clinical environments of all sizes. Both systems feature an enhanced design with smart automation and Konica Minolta's signature tube-mounted control interface. The KDR FMT is ideal for low-to mid-volume imaging needs, offering precision and ease of use. The KDR FMT+ adds automatic detector tracking, making it the perfect solution for mid-to high-volume practices seeking to boost throughput and enhance patient care. With advanced imaging software, flexible configurations, and immediate results, the KDR FMT and FMT+ empower clinicians to make **better decisions, sooner.**

IMAGING FLEXIBILITY, IMAGE RESOLUTION AND IMMEDIATE RESULTS.

The KDR FMT and FMT+ X-ray Systems.

TUBE-MOUNTED
CONTROL AND IMAGE DISPLAY


**MANUAL (FMT) or
AUTO (FMT+)**
TRACKING

**MULTIPLE TUBE/
GENERATOR**
COMBINATIONS

AeroDR®
14" X 17" AND 17" X 17"
FLAT PANEL DETECTOR OPTIONS

ELEVATING TABLE
700 LBS. CAPACITY

**FLOOR-MOUNTED
TUBE STAND**



EVERYTHING YOU EXPECT AND MORE.

SUPERIOR IMAGE QUALITY

The FMT Digital Radiography Systems captures high-resolution images in seconds. Combined with ULTRA imaging software, the result is a detailed visualization of both bone and soft tissue for a comprehensive diagnostic picture from a single exposure.

FLOOR-MOUNTED DESIGN

The Table and Wall Stand are designed to provide flexibility in a digital system that meets your clinical and economic needs.

WORKFLOW FLEXIBILITY

The Table and Wall Stand combination enables all imaging views required for general radiology, accommodating patients who are standing, sitting or lying on the table.



Specifications

Floor-Mounted Tube Stand	Collimator	Elevating Table
<ul style="list-style-type: none">FMT+ Auto Tracking to table or wall standFully counterbalanced tube armMechanical mid point, 90 and 180 degree detent	<ul style="list-style-type: none">Manual standard on FMT & FMT+Automatic collimator upgrade on FMT+ is available	<ul style="list-style-type: none">Motorized
<ul style="list-style-type: none">14" x 17" and 17" x 17" detector configurationsOptional rotating bucky for 14" x 17" detector	<ul style="list-style-type: none">Generator control from console and touchscreen	<ul style="list-style-type: none">4-way floating top700 lbs. lift capacity
<ul style="list-style-type: none">Longitudinal, transverse, vertical and rotational electromagnetic lock buttons	<ul style="list-style-type: none">Field lightCenter alignment light	<ul style="list-style-type: none">Foot controlsOptional table side control

BLUE MOON LIFECYCLE PLANS FOR X-RAY SYSTEMS. DO MORE. WORRY LESS.

Konica Minolta has developed a range of value-rich service solutions that prioritize helping our customers maximize their investments every day and respond immediately, or even proactively, to critical issues. Minimize downtime, maximize productivity, and mitigate risk with service solutions tailored to your needs.



**To learn more about
KDR FMT and KDR FMT+, contact your
Konica Minolta Sales Representative.**



KDR is a registered trademark of Konica Minolta, Inc. © 2026 Konica Minolta Healthcare Americas, Inc.



Konica Minolta Healthcare Americas, Inc.
411 Newark Pompton Turnpike
Wayne, New Jersey 07470
Tel: +1 (973) 633-1500
km.marketing@konicaminolta.com
healthcare.konicaminolta.us

M2405 1025 RevA