



# DIGITAL X-RAY PODIATRY SOLUTIONS

Revolutionary Imaging Systems



[2020imaging.net/podiatry](http://2020imaging.net/podiatry)



# Wireless and Tethered DR



## MOMENTUM DR 10X12 WIRELESS PANEL



MFG BY KONICA MINOLTA

- Wireless Interface
- Featherlight 3.30 lbs
- Cesium Scintillator
- 100 micron | 5 lp/mm
- Lithium Ion Capacitor  
(7-yr warranty)
- Cassette Sized
- 16-bit Images
- Rapid Charge

## UNIVERSAL DR 10X13 WIRELESS PANEL



UNIVERSAL DR

- Wireless Interface
- Lightweight 5 lbs
- Cesium Scintillator
- 100 micron | 5 lp/mm
- Lithium Ion Battery  
Extended Battery Life
- Active Image Area: 10"x13"
- 16-bit Images

## UNIVERSAL DR 1012V TETHERED PANEL



UNIVERSAL DR

- Tethered Interface
- Lightweight 5.5 lbs
- Cesium Scintillator
- 125 micron | 4 lp/mm
- Detachable & Durable Cable
- Cassette Sized
- 16-bit Images

## UNIVERSAL DR 14X17 WIRELESS PANEL



UNIVERSAL DR

- Wireless Interface
- Lightweight 7.3 lbs
- Cesium Scintillator
- 150 micron | 3.4 lp/mm
- Lithium Ion Battery
- Cassette Sized
- 16-bit Images

*Optional high resolution panel available*



# New and Retrofit X-Ray




## NEW PXS710D X-RAY

- High Frequency Output\*
- Dose Efficient/Ext Tube Life
- Superb Image Quality
- Bi-directional Tube Movement
- Integrated Cable Management
- Enhanced Patient Support
- UL Certified 450lbs weight limit
- Qualifies for ADA Tax Credit
- UL & CUL Listed
- Health Canada Approved

## RETROFIT ANY X-RAY

Wireless & Tethered DR Compatible  
All options qualify for ADA Tax Credit



Example of a converted X-Cel X-Ray, retrofitted to the BD-715 base/arm assembly with new 14x17 Universal DR panel.

**NOTE:** Existing Tube head & control utilized.



# Symmetry-Sync Software

## ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Background Submission
- Preset X-ray Exams
- Adult & Pediatric Settings
- Anatomical Program (APR)
- Fully Customizable
- Advanced Image Processing
- Image Enhancement Tools
- Touchscreen Enabled
- User Friendly Interface


## ADVANCED CLIENT VIEWER



### Multi-funtion Studylist

### Advanced Viewer

- EHR Integration Ready\*
- Podiatric Toolset Included
- Burn/Import CDs
- Fast, User Friendly Interface
- Custom Search/Query
- Image Fine Tuning
- Customized User Settings
- Exam Room Viewing



## ZERO FOOTPRINT WEB VIEWER

iPad®, iPhone®, Android™, tvOS, any mobile or smart device!

- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search
- Diagnostic DICOM Viewing







Advanced  
Superior  
Efficient  
Reliable  
Affordable

[2020imaging.net/podiatry](http://2020imaging.net/podiatry)

