



KONICA MINOLTA

**Contact:**

Dennis Puccio  
Konica Minolta Medical Imaging  
973-633-1500  
Web Address: [konicaminolta.com/medicalusa](http://konicaminolta.com/medicalusa)

**FOR IMMEDIATE RELEASE**

## **Konica Minolta Introduces Nano P+ Windows 7 Based CR System for Private Practices**

Wayne, NJ, August 7, 2014 – Konica Minolta Medical Imaging announced today the introduction of the Nano P+ CR system, a simple and easy-to-use solution designed for the private practice that provides the flexibility to upgrade to Digital Radiography. The Nano P+ offers enterprise level capabilities in a secure, Windows 7 based platform that enhances workflow and improves productivity.

Konica Minolta continues to provide customer focused Primary Imaging Solutions that offer clinical confidence and economic value with the introduction of the new Nano P+. The Nano P+ is a high performance, single bay CR system that is ideal for busy private practices. It provides an intuitive user interface with the CS-7 software and includes several menu and flow improvements. A more robust RIS/PACS interface includes modality work list functionality that enables the user to more easily customize exam tags or set-up multi-view tags, allowing for finer adjustments by the user.

With the Nano P+, practices can easily upgrade from CR to DR while maximizing value. The workstation can control more than one Konica Minolta digital imaging system, as well as the generator and collimator. With a compact design and flexible placement options, it is an ideal solution for private practices where space may be limited.

“With the launch of the Nano P+ CR system, Konica Minolta reaches new and current customers with solutions tailored to specific healthcare environments, and that are easily upgradeable for high economic value,” says Diane Hunt, Vice President of Marketing, Konica Minolta Medical Imaging. “This Windows 7 platform enhancement can help satisfy HIPAA requirements. Plus, with the launch of the Nano P+ with the CS-7NP

platform, Konica Minolta also continues to support its customer base by completing our transition from XP to Windows 7 systems.”

## About Konica Minolta Medical Imaging

Konica Minolta Medical Imaging is a world class provider and market leader in medical diagnostic primary imaging. With over 75 years of endless innovation, Konica Minolta is globally recognized as a leader providing cutting-edge technologies and comprehensive support aimed at providing real solutions to meet customer's needs. Konica Minolta Medical Imaging, headquartered in Wayne, NJ, is a unit of Konica Minolta, Inc. (TSE: 4902). For more information on Konica Minolta Primary Imaging Solutions, please visit [www.konicaminolta.com/medicalusa](http://www.konicaminolta.com/medicalusa).

<b>Company name</b>	KONICA MINOLTA, INC.
<b>Headquarters</b>	JP TOWER, 2-7-2 Marunouchi, Chiyoda-ku, Tokyo, Japan
<b>Founded</b>	December 1936
<b>FY 2013 Revenue</b>	\$9.4 Billion
<b>Number of employees</b>	Approx. 40,400 (2013)
<b>Business Lines</b>	The Konica Minolta Group operates in sectors ranging from business technologies, where our products are typified by MFPs (multi-functional peripherals), and Industrial Business (former Optics Business), where our products include pickup lenses for optical disks, and TAC film, a key material used in LCD panels, to healthcare, where we make digital X-ray diagnostic imaging systems.

### Konica Minolta Medical Imaging

411 Newark Pompton Turnpike; Wayne, NJ 07470 1-800-934-1034