

Real-time information for better decisions, sooner.



Get the secure, real-time information you need to maximize system utilization, minimize service interruptions, and evaluate personnel performance with the AeroRemote™ monitoring service for Konica Minolta AeroDR® Systems. This cloud-based, active monitoring solution detects and reports critical and ongoing performance issues, and provides automatic reports that you can review anywhere. This allows you to act on urgent situations immediately, and address potential system or user issues before they become problems.

Active Monitoring

Real-time monitoring of system parameters and events—from drop detection to component performance—lets you make better-informed support decisions, faster, while encouraging proper user handling. Better for your patients, better for your practice.

Live Notifications

Receive instant email alerts when urgent events, such as dropped panels, happen¹. Respond immediately and effectively, minimizing downtime and disruption.

Utilization Reporting

Monthly reports detailing system usage and personnel productivity for at-a-glance information helps you run a more efficient operation, support user best practices, and make the most of your AeroDR imaging assets. The reports provide details by system and by user:

- Panel drop counts
- Component health indicators
- OS notifications
- Exam counts
- Reject/repeat analysis
- Exposure level statistics

AeroRemote™: A complete active monitoring solution for your AeroDR® System.

Automated Panel Drop Detection, Diagnosis, and Replacement


With AeroRemote, you'll receive a real-time email notification when an AeroDR panel is dropped or bumped. Users are prompted by the AeroDR Control Station to run a panel check tool that gives an immediate result.

If the panel is damaged, a replacement is delivered via same-day or overnight shipping². Even installation is fast and easy. Simply insert the replacement AeroDR panel into the docking station, select the Quick Replacement option, and follow the instructions. In a matter of minutes you'll be ready for your next patient.



User-run Diagnostic Tool

Simplifies diagnosis with a user-driven panel check tool




Monthly Reports Accessible from Your Secure Portal

From panel drops and exposure levels to reject counts and more, a full menu of reports keeps you on top of system and operator productivity.



Rapid Replacement Service

Shortens the panel replacement time



1. Automated panel drop detection is not available on all models.
2. Panel replacement copays or other charges apply, depending on service contract level.

- Aero Remote requires system with CS-7 Version 1.21 or higher.

AeroDR is a registered trademark and AeroRemote is a trademark of Konica Minolta, Inc.
© 2016 Konica Minolta Healthcare Americas, Inc.



Konica Minolta Healthcare Americas, Inc.
411 Newark Pompton Turnpike
Wayne, New Jersey 07470
Tel: (973) 633-1500 Fax: (973) 523-7408
konicaminolta.com/medicalusa

M1007 1116 RevA