

Superior Imaging and Needle Visualization.



SONIMAGE® HS1 Compact Ultrasound System Transducers

Superior imaging with unique transducer technology, enhances clinical confidence.

The SONIMAGE® HS1 L18-4 is a broad frequency linear probe that offers both high sensitivity and deep penetration for scanning both superficial and deep joint structures.

The HL 18-4 hockey stick linear transducer is especially suited to MSK and pain management applications, reaching difficult to access areas with its small footprint and maneuverability.

SONIMAGE® HS1 Compatible Transducers

Konica Minolta transducers balance the trade-off between resolution and penetration, allowing superior imaging in both the near and far fields.



HL18-4 Linear Array

Application: MSK, Anesthesia, Pain
Frequency Range: 4 MHz to 18 MHz
Center Frequency: 10 MHz
Field of View: 25.6mm



L18-4 Linear Array

Application: MSK, Nerve, Pain Management
Frequency Range: 4 MHz to 18 MHz
Center Frequency: 10 MHz
Field of View: 38mm



L11-3 Linear Array

Application: MSK, Anesthesia, Pain
Frequency Range: 3 MHz to 11 MHz
Center Frequency: 7 MHz
Field of View: 38.4mm



L14-4 Linear Array

Application: Superficial, MSK, Vascular
Frequency Range: 4 MHz to 14 MHz
Center Frequency: 9 MHz
Field of View: 38.4mm



EC9-3 Endocavity Array

Application: Gynecology, Emergency
Frequency Range: 3 MHz to 9 MHz
Center Frequency: 5.7 MHz
Field of View: 176°



MC10-3 Micro Convex Array

Application: MSK, Anesthesia, Pain
Frequency Range: 3 MHz to 10 MHz
Center Frequency: 5.7 MHz
Field of View: 108°



S4-2 Phased Array

Application: Cardiac, Abdominal
Frequency Range: 2 MHz to 4 MHz
Center Frequency: 2.7 MHz
Field of View: 90°



C5-2 Convex Array

Application: Nerve
Frequency Range: 2 MHz to 5 MHz
Center Frequency: 3 MHz
Field of View: 51mm

Accessories

HL18-4: CIVCO® Latex-Free CIV-Flex™ Endocavity Probe Covers (#610-006)	MC10-3: CIVCO Latex-Free CIV-Flex™ Endocavity Probe Covers (#610-006)
L18-4: CIVCO Latex-Free CIV-Flex™ Endocavity Probe Covers (#610-542)	S4-2: CIVCO Latex-Free CIV-Flex™ Endocavity Probe Covers (#610-637)
L11-3: CIVCO Ultra-Pro II™ needle guides (#693-006), CIV-Flex™ Covers (#610-608)	C5-2: CIVCO Latex-Free CIV-Flex™ Endocavity Probe Covers (#610-637)
EC9-3: CIVCO Latex-Free CIV-Flex™ Endocavity Probe Covers (#610-006)	

© 2018 Konica Minolta Healthcare Americas, Inc. *Civco is a registered trademark of Civco Medical Solutions. Civ-Flex and Ultra-Pro II are trademarks of Civco Medical Solutions.



KONICA MINOLTA

Konica Minolta Healthcare Americas, Inc.
 411 Newark Pompton Turnpike
 Wayne, New Jersey 07470
 Tel: (973) 633-1500 Fax: (973) 523-7408
konicaminolta.com/medicalusa

M1025 0718 Rev B