

# Adjustable height, mobile cart for the J5 Ultrasound System.



## Portable Ergonomic Design

The J5 System can easily be mounted on this ergonomically designed, adjustable mobile cart.

Whether sitting or standing, the J5 mobile cart delivers a flexible, height-adjustable point of care cart. Even though it has a small footprint to fit tight spaces, this cart holds up to 23 lbs. A metal basket is included for gel, gloves, and other supplies.

- Locking mechanism secures the J5 System onto the cart
- Easily maneuvers via large 4" wide twin-wheel casters to the patient bedside, even over tricky thresholds.
- Ideal for point-of-care ultrasound



# Portable Ultrasound Cart for Point-of-Care Applications at the Bedside

With a smart, compact and height adjustable design, the **Konica Minolta J5 Mobile Ultrasound Cart** is built with MSK physicians' requirements in mind. Carry all your ultrasound supplies right to your patient on this compact, easy-to-steer pole cart.

## Cart Specifications:

- Weight: (Without Device): 10 lb
- Load range: 5 to 23 lb/2.3 to 10.4 kg
- Height adjustment range: 11"/28 cm
- Minimum height: 34.5"/87.6 cm
- Maximum height: 45.5"/115.6 cm
- 21"/53.5 cm diameter aluminum base
- 10 lb/4.5 kg counterweight

## Cart Features:

- 4" / 10.2 cm diameter dual wheel nylon locking casters for easy transition over uneven surfaces
- Gas spring assisted for easy height adjustment
- Foot pedal or hand knob for height adjustment

## Included with Cart:

- All required hardware
- Directions for use (DU)
- Non-standard tools required for assembly

## Order Information:

- Adjustable height, mobile cart for J5 System (# 5905400B)



*Adjustable height, mobile cart (#5905400B) with J5 System (sold separately)*



**KONICA MINOLTA**

**Konica Minolta Healthcare Americas, Inc.**

411 Newark Pompton Turnpike

Wayne, New Jersey 07470

Tel: (973) 633-1500 Fax: (973) 523-7408

[konicaminolta.com/medicalusa](http://konicaminolta.com/medicalusa)

M1108 517 Rev A