

CS-7 Universal Controller



CS-7 Universal Controller is an intuitive PC with a multi-touch flat screen display designed to provide maximum performance with ease of use and clinical confidence in the radiology department. Its sophisticated multiple functions will provide improved workflow of X-ray examinations. Obtain patient/examination information from RIS with Modality Worklist, select a patient, perform the exams, verify images, and send to PACS for radiologist's review/diagnosis. CS-7 can control AeroDR® detectors and Xpress CR® readers and provide a path to selected X-ray controllers. Konica Minolta's *Hybrid (H) Processing* is included with every CS-7 Controller.

Hybrid (H) Processing is a unique image processing technique developed to provide superior image quality and improved patient diagnosis. This enables the technologist to acquire, view, and complete a study in a smooth, time saving workflow, spending less time at the image Control Station and more time with patient care, resulting in improved workflow.

Quick Preview and smart User Interface

After exposure, a preview image immediately appears on the CS-7 console in *less than two seconds*. The powerful user-friendly interface provides proprietary functions aiding in quick diagnosis for better patient care.

