



Motion / Lung Detection

Motion Blur




Automatic Warning

- Possible motion blur detected
- Circled with **Orange frame**.
- Magnifier button displayed.




Anatomy Cutoff



Automatic Warning

- Anatomy cutoff detected
- Circled with **Red frame**.



Optional/additional CS-7 Software application

Motion/Lung Detection is software that automatically highlights potential motion artifacts and possible anatomy cut-off on a chest X-ray, alerting the technologist to potential issues with the acquisition. The software further enhances clinical confidence and simplifies workflow with a magnifier button for immediate assessment, which maximizes reliability.