

Study Combine / Move



Optional/additional CS-7 Software application

Study Combine/Move provides improved patient experience. It provides the technologist the ability to move an existing patient study (Study A), which has been originally registered under a particular patient name, to a different patient name received via modality worklist (Study B), while still retaining the registration information from Study A. This feature is also a useful tool when correcting a patient study that has been incorrectly mislabeled, or in trauma environments where the patient identity is not known at the time of the initial image acquisition.

Examples:

- Orders for the same patient/body part are issued from more than one doctor on the same day
- If images are mixed between patients, the technologist has the ability to move the study and change patient information
- Images taken with a provisional ID# can be moved to a regular work request number
- Images can be moved from one study to another

Example 1: Using a temporary order, an exam can be performed on an unknown patient (e.g. emergency)

- **Benefit** – Update later, with correct, confirmed patient information

Example 2: Performing an exam on “Patient A,” but using the order of “Patient B” by mistake

- **Benefit** - Provides ability to update/correct patient information