


A radiography solution, on the move.



State-of-the-art imaging software supports detailed bone and soft tissue visualization

The mKDR[®] Digital Radiography System brings advanced digital X-ray capabilities anywhere at your facility.

Clinicians will have the flexibility, image resolution and immediate results needed to make informed decisions faster, boosting productivity and increasing patient satisfaction.

- Superior image quality
- Compact footprint
- Versatile design

A system that goes anywhere you go.

Superior image quality

The mKDR Essential includes a 14" x 17" detector that captures digital images in seconds. When combined with Ultra imaging software modified for mobile use, detailed bone and soft-tissue images are provided from a single exposure.

Space-saving design

The mKDR features a collapsing column and swivel arm to maximize visibility and positioning without sacrificing clinical utility.

Engineered for versatility

A full range of motion enables imaging at bedside, NICU, ICU, ER and even as a backup system for exam rooms.

Designed for today's environment

AeroDR® Digital Radiography Solutions combine performance and productivity to help your imaging practice run more smoothly and efficiently.

- 14"x17"
- Multiple configurations
- Fully integrated

1. Collapsible Column

2. Convenient Swivel Arm Provides A Wide Range of Motion

3. Power Driven

4. In-Bin Detector Charging



FEATURES/SPECIFICATIONS		
Movement	Collimator	Options
Power drive	Manual Control	Multiple Tube/Generator Combinations
Manual Rotation with Detents: 30°~120° - Detents on 30,60,90 and 120 degrees	Accessory Rails (Cones, Filters, DAP)	Pediatric Imaging Package
Collapsible Column	Field Light	
Tube Rotation ±90°	Six Blades	

Value, confidence and worry-free operation.

Konica Minolta's high customer satisfaction rating is based on supporting you throughout the entire lifecycle of your system, not just when it breaks. Blue Moon's wider range of benefits delivers economical value and clinical confidence.

- Minimize downtime
- Maximize productivity
- Eliminate Risk



BLUE MOON

Select the right plan to meet your needs.

**Blue Moon
Elite**

**Blue Moon
Select**

**Blue Moon
Basic**

©2020 Konica Minolta Healthcare Americas, Inc.



KONICA MINOLTA

Konica Minolta Healthcare Americas, Inc.

411 Newark Pompton Turnpike
Wayne, New Jersey 07470
Tel: (973) 633-1500 Fax: (973) 523-7408
konicaminolta.com/medicalusa

M1282 0218 RevB